

執行 JpegUploader



輸入 HIS 帳號登入

A dialog box titled '登入陽明HIS帳號密碼' (Login阳明HIS account password). It contains two text input fields: '使用者名稱' (Username) and '使用者密碼' (Password). Below the fields is a checkbox labeled '記住帳號密碼' (Remember account password). At the bottom are two buttons: '登入(L)' (Login) and '取消(C)' (Cancel).

輸入搜尋條件

選好日期範圍，區間內的資料都會顯示

A search criteria form with the title '1. 查詢檢查單' (1. Query check sheet). It includes a date range selector set to '2023/03/21' to '2023/04/21', a '類別' (Category) dropdown, and a 'Modality' dropdown. There are also fields for '病歷號' (Medical record number), '申請單號' (Application number), and a '是否已上傳' (Whether uploaded) dropdown set to 'Y'. A '查詢' (Query) button is on the right.

輸入病歷號前面補 0 到 10 碼

The same search criteria form as above, but the '病歷號' (Medical record number) input field is highlighted with a red box, indicating where to enter the ID.

輸入檢驗單號

The same search criteria form as above, but the '申請單號' (Application number) input field is highlighted with a red box, indicating where to enter the application ID.

是否以上傳，Y 是會顯示已上傳的資料，N 是顯示還未上傳資料

(1936,1056)

檔案(F) 設定(S)

1.查詢檢查單

日期範圍: 2023/03/21 ~ 2023/04/21 類別: Modality:

病歷號: 申請單號: 是否已上傳: Y 查詢

影像上傳程式(UPN)

1.查詢檢查單

日期範圍: 2023/02/28 ~ 2023/04/24 類別: Modality:

病歷號: 申請單號: 是否已上傳: N 查詢


申請單號	CODE	ACC NO	病歷號碼	病患姓名	Modality	OrderDate
B12120304001	心電圖	B12120304001	0000082019	洪耀輝	SC	1120304
B13120317001	其他超音波	B13120317001	0000165260	鄭培權	XXXXXXXX	1120317
B13120306001	其他超音波	B13120306001	0000143446	陳威達	XXXXXXXX	1120306
B13120306002	其他超音波	B13120306002	0000163622	陳宇洋	XXXXXXXX	1120306
B13120308002	其他超音波	B13120308002	0000022706	邱華燦	XXXXXXXX	1120308
B13120308003	其他超音波	B13120308003	0000165080	梁秀梅	XXXXXXXX	1120308
B13120308004	其他超音波	B13120308004	0000043292	林素馨	XXXXXXXX	1120308

2.選擇影像檔案

④ JPG ○ PDF

3. 影像縮圖

4.上傳影像 開始上傳(U)



1. 點選要上傳的單號，黃色表示已選取此筆資料
2. 點選+號選取路徑裡的圖片
3. 要選擇多張圖片時按住鍵盤 Ctrl+滑鼠點選圖片檔案

申請單號	CODE	ACC NO
B12120304001	心電圖	B12120304001
B13120317001	其他超音波	B13120317001
B13120306001	其他超音波	B13120306001
B13120306002	其他超音波	B13120306002
B13120308002	其他超音波	B13120308002
B13120308003	其他超音波	B13120308003
B13120308004	其他超音波	B13120308004

2.選擇影像檔案

JPG  PDF  

序號	檔名	檔案位置
1	02.jpg	C:\Users\MIS\Desktop\EKG上傳步驟
2	04.jpg	C:\Users\MIS\Desktop\EKG上傳步驟
3	05.jpg	C:\Users\MIS\Desktop\EKG上傳步驟

4.上傳影像 

1. 選取的圖片檔案會顯示出來
2. 點選開始上傳